

Probe Wipes

9150 & 9152

Probe Wipes help meet HACCP guidelines and are an ideal way of cleaning thermometer probes between temperature checks to avoid cross-contamination. The 9152 tub and the 9150 box of individual packets each contain 200 wipes, each wipe consists of 70% Isopropyl alcohol. Sanitizing probes before and after use is recommended.

9152 Probe Wipe Tub



Specifications

9150 Probe Wipe Box

- 200 Count
- Ingredients: 70% Isopropyl alcohol
- Wipe Dimensions: 2" x 2" (50 mm x 50 mm)
- One year warranty

9152 Probe Wipe Tub

- 200 Count
- Ingredients: 70% Isopropyl alcohol
- Wipe Dimensions: 3" x 5" (76 mm x 127 mm)
- One year warranty

Packaging

9150 Probe Wipe Box

- Individual Box Package Weight: 7 oz (198 g)
- Boxes Per Carton: 10
- Carton Dimensions: 11.5" x 8.88" x 6.25" (292 mm x 225 mm x 159 mm)
- Package Cube: 0.37

9152 Probe Wipe Tub

- Individual Tub Package Weight: 6 oz (170 g)
- Tubs Per Carton: 24
- Carton Dimensions: 12.25" x 9.5" x 9.75" (311 mm x 241 mm x 248 mm)
- Package Cube: 0.65

9150 Probe Wipe Box



Accuracy to the Highest Degree

Cooper-Atkins Corporation • 800-835-5011 • 860-347-2256 • www.cooper-atkins.com

67-1307
V0213